

Product Datasheet

SARS-CoV-2 (2019-nCoV) S RBD-SD1 protein (orb669698)

Catalog Number	orb669698
Category	Proteins
Description	Recombinant SARS-CoV-2 (2019-nCoV) S Protein RBD-SD1 is produced by our Mammalian expression system and the target gene encoding Arg319-Ser591 is expressed with a 6His tag at the C-terminus.
Tag	C-His
UniProt ID	QHD43416.1
Expression System	Human Cells
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Biological Origin	SARS-CoV-2
Expression Region	Arg319-Ser591
Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Dry Ice Shipping	Please note: This product requires shipment on dry ice. A dry ice surcharge will apply.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)