

Product Datasheet

Goat Anti-Rabbit IgG (H+L) HRP (orb670232)

Catalog Number	orb670232
Category	Antibodies
Description	Peroxidase-conjugated AffiniPure Goat Anti-rabbit IgG (H+L)
Species/Host	Goat
Conjugation	HRP
Reactivity	Rabbit
Concentration	1 mg/ml
Buffer/Preservatives	0.02% thimerosal. 1mg/ml in PBS with 0.02% thimerosal, 50% glycerol, pH 7.4
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Immunogen	Rabbit IgG (H+L)
Tested applications	ELISA, WB
Dilution range	EIA and Western blots: 1:5,000-1:100,000. Optimal working dilutions must be determined by end user.
Antibody Type	Secondary Antibody
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)