

Product Datasheet

Mouse Anti-Goat IgG (H+L) Rhodamine (TRITC) (orb670242)

Catalog Number	orb670242
Description	Fluorophore: Tetramethyl rhodamine 5(and6)-isothiocyanate(TRITC) Amax=550nm; Emax=570nm Fluorophore/ProteinAbsorbanceRatio: A550/A280=0.48
Species/Host	Mouse
Reactivity	Goat
Conjugation	Rhodamine
Tested Applications	FC, IF
Immunogen	Goat IgG(H+L)
Preservatives	0.05% Sodium Azide. 0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6,containing 1% rAlbumin,50% glycerol
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at - 20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Antibody Type	Secondary Antibody
Dilution Range	1:50 -1:200 for most applications
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)