

Product Datasheet

Goat Anti-Rabbit IgG (H+L) Biotin (orb670245)

Catalog Number	orb670245
Category	Antibodies
Description	Biotin-SP-conjugated AffiniPure Goat Anti-Rabbit IgG (H+L)
Species/Host	Goat
Conjugation	Biotin
Reactivity	Rabbit
Concentration	1 mg/ml
Buffer/Preservatives	0.05% Sodium Azide. 0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, containing 1% rAlbumin, 50% glycerol
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Immunogen	Rabbit IgG (H+L)
Tested applications	ELISA, FC, ICC, IHC, WB
Dilution range	1:20000-1:400,000 for ELISA and Western blotting using enzyme-conjugated streptavidin; 1:500-1:5,000 for enzyme immunohisto/cytochemistry; 1:200-1:1,000 for flow cytometry and fluorescence immunohisto/cytochemistry
Antibody Type	Secondary Antibody
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Expiration Date 12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)