

## Product Datasheet

### Rabbit Anti-Goat IgG(H+L) Cy3 (orb670257)

<b>Catalog Number</b>	orb670257
<b>Category</b>	Antibodies
<b>Description</b>	Fluorophore: Cy3 bis-NHS ester Amax=550 nm; Emax=570 nm Fluorophore/Protein: A552/A280=2.3
<b>Species/Host</b>	Rabbit
<b>Conjugation</b>	Cy3
<b>Reactivity</b>	Goat
<b>Concentration</b>	1 mg/ml
<b>Buffer/Preservatives</b>	0.05% Sodium Azide. 0.01M Sodium phosphate, 0.25M NaCl, 50% Glycerol, pH 7.6
<b>Purification</b>	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
<b>Immunogen</b>	Goat IgG(H+L)
<b>Tested applications</b>	FC, IF
<b>Dilution range</b>	1:100 - 1:800 for most applications
<b>Antibody Type</b>	Secondary Antibody
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Expiration Date</b>	12 months from date of receipt.

**Biorbyt Ltd.**

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

**Biorbyt LLC.**

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)