

Product Datasheet

HAKAI (C) Antibody, Rabbit Polyclonal (orb67187)

Catalog Number	orb67187
Category	Antibodies
Description	Rabbit polyclonal to HAKAI.HAKAI is also known as Casitas B-lineage lymphoma-transforming sequence-like protein 1 (CBLL1),which is an E3 ubiquitin-protein ligase for the E-cadherin (CDH1) complex and mediates itsubiquitination, endocytosis, and degradation in the lysosomes. CDH1 is endocytosed as aconsequence of tyrosine phosphorylation and ubiquitination. HAKAI interacts with tyrosinephosphorylatedCDH1 and mediates this process. HAKAI contains a RING-finger domain and is alsothought to have a role in the control of cell proliferation.
Clonality	Polyclonal
Species/Host	Rabbit
Conjugation	Unconjugated
Reactivity	Human
Buffer/Preservatives	0.01% NaN3
Purity	Affinity purification
Immunogen	C-terminal region of human HAKAI
UniProt ID	Q75N03
MW	54.5
Tested applications	WB
Dilution range	WB:1:500-1:1000

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)


Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.



Western blot analysis of total cell extracts from human HeLa using HAKAI antibody

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)