

## Product Datasheet

### ICAM1 Antibody (orb675641)

|                             |                                                                                                                                                                                                                                                                                                                                                                               |
|-----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog Number</b>       | orb675641                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Category</b>             | Antibodies                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Description</b>          | This ICAM1 Antibody is an unconjugated polyclonal product. It targets ICAM1 using a synthesized peptide derived from Human ICAM-1 around the non-phosphorylation site of Y512 as the immunogen. This antibody is suitable for ELISA, WB. Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.. |
| <b>Target</b>               | ICAM1                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Clonality</b>            | Polyclonal                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Species/Host</b>         | Rabbit                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Isotype</b>              | IgG                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Conjugation</b>          | Unconjugated                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Reactivity</b>           | Human                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Form/Appearance</b>      | Liquid                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Buffer/Preservatives</b> | Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.                                                                                                                                                                                                                                                                                                  |
| <b>Purification</b>         | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.                                                                                                                                                                                                                                                         |
| <b>Immunogen</b>            | Synthesized peptide derived from Human ICAM-1 around the non-phosphorylation site of Y512.                                                                                                                                                                                                                                                                                    |
| <b>UniProt ID</b>           | <b>P05362</b>                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Tested applications</b>  | ELISA, WB                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Dilution range</b>       | WB: 1:500-1:2000, ELISA: 1:10000                                                                                                                                                                                                                                                                                                                                              |

#### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

#### Biorbyt LLC.


68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

**Storage**


Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

**Note**

For research use only



Western blot analysis of PC3 22RV1 HELA cells using ICAM1 antibody



Western blot analysis of HepG2 cells using ICAM1 antibody

**Biorbyt Ltd.**

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC.**

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)