

Product Datasheet

Goat Anti-Mouse IgG(H+L) Antibody;HRP conjugated (orb688922)

Catalog Number	orb688922
Category	Antibodies
Description	This Goat Anti-Mouse IgG(H+L) Antibody;HRP conjugated is an HRP polyclonal product, using a Mouse IgG(H+L) as the immunogen. This antibody is suitable for ELISA, FC, IF, IHC, IP, WB. Purification: The antibody was affinity-purified from goat antiserum by affinity-chromatography using specific immunogen..
Clonality	Polyclonal
Species/Host	Goat
Isotype	IgG
Conjugation	HRP
Reactivity	Mouse
Form/Appearance	Liquid
Buffer/Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Purification	The antibody was affinity-purified from goat antiserum by affinity-chromatography using specific immunogen.
Immunogen	Mouse IgG(H+L)
Tested applications	ELISA, FC, IF, IHC, IP, WB
Dilution range	WB:1:5000-20000, IHC-P:1:200-2000, IF:1:200-2000, IP:1:200-1:2000, FCM:1:200-2000, ELISA: 1:5000-20000
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)