

Product Datasheet

Goat Mouse IgG(H+L) Antibody, FITC conjugated (orb688924)

Catalog Number	orb688924
Category	Antibodies
Description	This Goat Anti-Mouse IgG(H+L) Antibody;FITC conjugated is an FITC polyclonal product, using a Mouse IgG(H+L) as the immunogen. This antibody is suitable for FC, IF, IHC. Purification: The antibody was affinity-purified from goat antiserum by affinity-chromatography using specific immunogen..
Clonality	Polyclonal
Species/Host	Goat
Isotype	IgG
Conjugation	FITC
Reactivity	Mouse
Form/Appearance	Liquid
Buffer/Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Purification	The antibody was affinity-purified from goat antiserum by affinity-chromatography using specific immunogen.
Immunogen	Mouse IgG(H+L)
Tested applications	FC, IF, IHC
Dilution range	WB:1:5000-20000, IHC-P:1:200-2000, IF:1:200-2000, IP:1:200-1:2000, FCM:1:200-2000, ELISA: 1:5000-20000

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)