

Product Datasheet

Goat Anti-Rabbit IgG(H+L) Antibody (orb688929)

Catalog Number	orb688929
Category	Antibodies
Description	This Goat Anti-Rabbit IgG(H+L) Antibody is an unconjugated polyclonal product, using a Rabbit IgG(H+L) as the immunogen. This antibody is suitable for ELISA, FC, IF, IHC, IP, WB. Purification: The antibody was affinity-purified from goat antiserum by affinity-chromatography using specific immunogen..
Clonality	Polyclonal
Species/Host	Goat
Isotype	IgG
Conjugation	Unconjugated
Reactivity	Rabbit
Form/Appearance	Liquid
Buffer/Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Purification	The antibody was affinity-purified from goat antiserum by affinity-chromatography using specific immunogen.
Immunogen	Rabbit IgG(H+L)
Tested applications	ELISA, FC, IF, IHC, IP, WB
Dilution range	WB:1:5000-20000, IHC-P:1:200-2000, IF:1:200-2000, IP:1:200-1:2000, FCM:1:200-2000, ELISA: 1:5000-20000
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)