

Product Datasheet

Affinipure Goat Rabbit IgG(H+L), HRP conjugated Antibody (orb688932)

Catalog Number	orb688932
Category	Antibodies
Description	This HRP-conjugated Affinipure Goat Anti-Rabbit IgG(H+L) is an HRP polyclonal product, using a Rabbit IgG(H+L) as the immunogen. This antibody is suitable for ELISA, WB. Purification: The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads..
Clonality	Polyclonal
Species/Host	Goat
Isotype	IgG
Conjugation	HRP
Reactivity	Rabbit
Form/Appearance	Liquid
Buffer/Preservatives	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6
Purification	The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Immunogen	Rabbit IgG(H+L)
Tested applications	ELISA, WB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Expiration Date

12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)