

## Product Datasheet

# NUP153 Recombinant Monoclonal Antibody (orb689152)

|                            |                                                                                                                                         |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog Number</b>      | orb689152                                                                                                                               |
| <b>Description</b>         | Rabbit monoclonal antibody to NUP153                                                                                                    |
| <b>Reactivity</b>          | Human                                                                                                                                   |
| <b>Conjugation</b>         | Unconjugated                                                                                                                            |
| <b>Tested Applications</b> | ELISA, IF, IHC                                                                                                                          |
| <b>Immunogen</b>           | A synthesized peptide derived from human Nup153                                                                                         |
| <b>Target</b>              | NUP153                                                                                                                                  |
| <b>Preservatives</b>       | Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.                                       |
| <b>Form/Appearance</b>     | Liquid                                                                                                                                  |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>                | For research use only                                                                                                                   |
| <b>Isotype</b>             | Rabbit IgG                                                                                                                              |
| <b>Clonality</b>           | Monoclonal                                                                                                                              |
| <b>Clone Number</b>        | 8D4                                                                                                                                     |
| <b>Antibody Type</b>       | Recombinant Antibody                                                                                                                    |
| <b>Uniprot ID</b>          | <b>P49790</b>                                                                                                                           |
| <b>Dilution Range</b>      | IHC-P:1:50-1:200, IF:1:20-1:200                                                                                                         |


### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+441223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)


### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

**Expiration Date** 12 months from date of receipt.



Immunofluorescence analysis of Hela cells using NUP153 antibody



Immunohistochemical staining of human breast cancer using NUP153 antibody

---

**Biorbyt Ltd.**

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

**Biorbyt LLC.**

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)