

Product Datasheet

Anti-SARS-CoV-2 (CB6 biosimilar) mAb (orb689316)

Catalog Number	orb689316
Description	Anti-SARS-CoV-2 (CB6 biosimilar) mAb
Species/Host	Human
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Target	SARS
Form/Appearance	Powder
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Isotype	IgG1
Clonality	Monoclonal
Antibody Type	Biosimilar Antibody
Uniprot ID	P0DTC2
Dilution Range	ELISA 1:5000-10000


Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Expiration Date 12 months from date of receipt.



ELISA plate pre-coated by 2 µg/ml (100 µL/well) Human S-RBD, mFc-His tagged protein orb689446 can bind Anti-SARS-CoV-2 Neutralizing antibody in a linear range of 0.24-15.62 ng/ml.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+441223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)