

Product Datasheet

BX430 (orb692024)

Description	BX430 is a potent and selective noncompetitive allosteric human P2X4 receptor channels antagonist with an IC ₅₀ of 0.54 μM. BX430 has species specificity. BX430 is used for chronic pain and cardiovascular disease.
Note	For research use only
Formula	C ₁₅ H ₁₅ Br ₂ N ₃ O
MW	413.11
CAS Number	[688309-70-8]
SMILES	O=C(NC1C(Br)=CC(C(C)C)=CC=1Br)NC1C=NC=CC=1
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)