

Product Datasheet

Donkey anti-Sheep IgG (H&L), Biotin conjugated (orb699951)

Catalog Number	orb699951
Description	Donkey polyclonal antibody against Donkey anti-Sheep IgG (H&L) conjugated to Biotin
Species/Host	Donkey
Reactivity	Sheep
Conjugation	Biotin
Tested Applications	ELISA, IHC, WB
Immunogen	Purified sheep IgG, whole molecule
Form/Appearance	Liquid
Storage	Store non-diluted antibody at 2-8°C. Shelf life is 1 year from date of receipt.
Note	For research use only
Application notes	No reactivity is observed to non-immunoglobulin mouse serum proteins based in immunoelectrophoresis. BSA and milk have to be replaced by other blocking reagents, like donkey serum or commercial formulations which are free from bovine IgG.
Clonality	Polyclonal
Purity	Immunogen purified donkey IgG.
Dilution Range	1 : 2000-1: 20 000 (ELISA), 1 : 500 -1 : 5000 (IHC), 1 : 2000-1: 20 000 (WB)
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)