

## Product Datasheet

### Goat anti-Rabbit IgG (H&L), F(ab)'2 fragment, TRITC conjugated (orb699973)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                  |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>       | Goat polyclonal antibody against Goat anti-Rabbit IgG (H&L) conjugated to TRITC                                                                                                                                                                                                                                                                                                                                                  |
| <b>Species/Host</b>      | Goat                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Conjugation</b>       | TRITC                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Form/Appearance</b>   | Liquid                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Storage</b>           | Store at 2-8°C. Shelf life of this product is one year from date of receipt.                                                                                                                                                                                                                                                                                                                                                     |
| <b>Note</b>              | For research use only                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Application notes</b> | This antibody reacts with the heavy chains on rabbit IgG and with the light chains on all rabbit immunoglobulins based on immunoelectrophoresis. Minimum cross-reactivity is observed to non-immunoglobulin rabbit serum proteins based on immunoelectrophoresis. Antibody is supplied in 10 mM sodium phosphate, 150 mM sodium chloride, pH 7.2, 1% (w/v) BSA, protease/IgG free and 0.05 % (w/v) sodium azide as preservative. |
| <b>Clonality</b>         | Polyclonal                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Purity</b>            | Immunogen affinity purified IgG, F(ab)'2 fragment.                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Dilution Range</b>    | The optimal working dilution should be determined by the investigator, (ICC), immunohistochemistry (IHC) (Flow cyt)                                                                                                                                                                                                                                                                                                              |
| <b>Expiration Date</b>   | 12 months from date of receipt.                                                                                                                                                                                                                                                                                                                                                                                                  |

#### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

#### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)