

Product Datasheet

Goat Rabbit IgG (H&L), F(ab)'2 fragment, TRITC conjugated, min, cross-reactivity to bovine,human, mouse IgG/serum Antibody (orb700002)

Catalog Number	orb700002
Category	Antibodies
Description	Goat polyclonal antibody against Goat anti-Rabbit IgG (H&L) conjugated to TRITC
Clonality	Polyclonal
Species/Host	Goat
Conjugation	TRITC
Reactivity	Rabbit
Form/Appearance	Liquid
Purification	Immunogen affinity purified IgG, F(ab)'2 fragment.
Dilution range	The optimal working dilution should be determined by the investigator.
Application notes	This antibody reacts with the heavy chains on rabbit IgG and with the light chains on all rabbit immunoglobulins based on immunoelectrophoresis.Minimum cross-reactivity is observed to non-immunoglobulin rabbit serum proteins and bovine,human,mouse IgG/serum based on immunoelectrophoresis.Antibody is supplied in 10 mM sodium phosphate, 150 mM sodium chloride, pH 7.2, 1% (w/v) BSA, protease/IgG free and 0.05% (w/v) sodium azide as preservative
Storage	Store at 2-8°C. Shelf life of this product is one year from date of receipt.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Expiration Date 12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)