

Product Datasheet

Goat Rabbit IgG (H&L) HRP conjugated, min,cross-reactivity to bovine/Human/mouse IgG/serum Antibody (orb700290)

Catalog Number	orb700290
Category	Antibodies
Description	Goat polyclonal antibody against Goat anti-Rabbit IgG (H&L) HRP conjugated conjugated to HRP
Clonality	Polyclonal
Species/Host	Goat
Conjugation	HRP
Reactivity	Rabbit
Form/Appearance	Lyophilized
Reconstitution	For reconstitution add 1,1 ml of sterile water, Let it stand 30 minutes at room temperature to dissolve, Prepare fresh working dilutions daily
Purification	Immunogen affinity purified goat IgG.
Immunogen	purified Rabbit IgG, whole molecule
Tested applications	ELISA, IHC, WB
Dilution range	1 : 500 -1 : 5000 (ELISA), 1 : 500 -1 : 5000 (IHC), 1 : 500 -1 : 5000 (WB)
Application notes	Based on immunoelectrophoresis no reactivity is observed to:non-immunoglobulin rabbit serum proteinsserum from bovine, human or mouseIgG from human or mouse

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Storage

Store lyophilized material at 2-8°C. For long time storage after reconstitution, dilute the antibody solution with glycerol to a final concentration of 50% glycerol and store as liquid at -20°C, to prevent loss of enzymatic activity. For example, if you have reconstituted 1 mg of antibody in 1,1 ml of sterile water add 1,1 ml of glycerol. Such solution will not freeze in -20°C. If you are using a 1:5000 dilution prior to diluting with glycerol, then you would need to use a 1:2500 dilution after adding glycerol. Prepare working dilution prior to use and then discard, Be sure to mix well but without foaming.

Note

For research use only

Expiration Date

12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)