

Product Datasheet

Chicken Rabbit IgG (H&L), Affinity purified, Unconjugated, min, cross-reactivity to human, mouse IgG Antibody (orb700355)

Catalog Number	orb700355
Category	Antibodies
Description	Gallus polyclonal antibody to Gallus anti-Rabbit IgG (H&L)
Clonality	Polyclonal
Species/Host	Gallus
Conjugation	Unconjugated
Reactivity	Rabbit
Form/Appearance	Liquid
Purity	Affinity purified chicken IgY
Purification	Immunogen affinity purified chicken IgY.
Immunogen	Purified Rabbit IgG
Tested applications	ELISA
Dilution range	The optimal working dilution should be determined by the investigator

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Application notes

Antibody reacts with the heavy chains on rabbit IgG and with the light chains on all rabbit immunoglobulins based on immunoelectrophoresis. No reactivity is observed to non-immunoglobulin rabbit serum and human or mouse serum or IgG based on immunoelectrophoresis. Antibody is supplied in 10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2, 0.05 % (w/v) of Sodium azide as preservative. Antibody concentration is 4.5 mg/ml

Storage

Store non-diluted antibody at 2-8°C, For storage at -20°C dilute antibody solution with an equal volume of glycerol to obtain final glycerol concentration of 50 % to prevent loss of enzymatic activity, Such solution will not freeze in -20°C, If you are using a 1:5000 dilution prior to diluting with glycerol, then you would need to use a 1:2500 dilution after adding glycerol, Prepare working dilution prior to use and then discard, Be sure to mix well but without foaming,!!!Store non-diluted antibody at 2-8°C, For storage at -20°C dilute with an equal volume of glycerol to prevent loss of enzymatic activity, Prepare working dilution prior to use and then discard, Shelf life is 1 year from date of receipt,

Note

For research use only

Expiration Date

12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558