

Product Datasheet

Goat Donkey IgG (H&L), Unconjugated Antibody (orb700464)

Catalog Number	orb700464
Category	Antibodies
Description	Goat polyclonal antibody to Goat anti-Donkey IgG (H&L)
Clonality	Polyclonal
Species/Host	Goat
Conjugation	Unconjugated
Reactivity	Donkey
Form/Appearance	Liquid
Purity	Affinity purified
Purification	Immunogen affinity purified goat IgG.
Immunogen	Purified Donkey IgG, whole molecule.
Tested applications	IA
Dilution range	The optimal working dilution should be determined by the investigator
Application notes	Antibody purity is > 95% based on SDS-PAGE. Affinity purified using solid phase Donkey IgG. Antibody is supplied in 10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2. 0.05% (w/v) Sodium Azide is added as a preservative. Based on IEP, this antibody reacts with: • heavy (γ) chains on donkey IgG • light chains on all donkey immunoglobulins. Based on IEP, no reactivity is observed to: • non-immunoglobulin donkey serum immunoglobulins

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Storage	Store at 2-8 °C,
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)