

Product Datasheet

Goat Human IgG (H&L), Biotin conjugated, min, cross-reactivity to rat IgG Antibody (orb700475)

Catalog Number	orb700475
Category	Antibodies
Description	Goat polyclonal antibody against Goat anti-Human IgG (H&L conjugated to Biotin
Clonality	Polyclonal
Species/Host	Goat
Conjugation	Biotin
Reactivity	Human
Form/Appearance	Liquid
Purification	Immunogen affinity purified goat IgG.
Immunogen	Purified Human IgG, whole molecule
Tested applications	FC, IL
Dilution range	1 : 20-1 : 2000 (Flow cyt), (IL)
Application notes	Antibody purity is > 95% based on SDS-PAGE.Affinity purified using solid phase Human IgG .Antibody is supplied in 10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2, 1% (w/v) BSA, Protease/IgG free. 0.05% (w/v) Sodium Azide is added as a preservative.Based on IEP, this antibody reacts with: • heavy (γ) chains on human IgG • light chains on all human immunoglobulins.Based on IEP, no reactivity is observed to: • non-immunoglobulin human serum immunoglobulins • rat IgG
Storage	Store at -20° C.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)