

## Product Datasheet

# Fbx32 Recombinant Rabbit Monoclonal Antibody (orb704278)

|                      |                                                                        |
|----------------------|------------------------------------------------------------------------|
| Catalog Number       | orb704278                                                              |
| Category             | Antibodies                                                             |
| Description          | Fbx32 Recombinant Rabbit Monoclonal Antibody                           |
| Target               | FBXO32                                                                 |
| Clonality            | Recombinant                                                            |
| Species/Host         | Rabbit                                                                 |
| Isotype              | IgG                                                                    |
| Conjugation          | Unconjugated                                                           |
| Reactivity           | Mouse, Rat                                                             |
| Predicted Reactivity | Human                                                                  |
| Form/Appearance      | Liquid                                                                 |
| Concentration        | 1mg/ml                                                                 |
| Buffer/Preservatives | 0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol. |
| Purification         | Affinity purified by Protein A                                         |
| Immunogen            | A synthesized peptide derived from human FBXO32 (20-84aa)              |
| UniProt ID           | <b>Q969P5</b>                                                          |
| MW                   | 42 kDa                                                                 |

### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom


Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

|                            |                                                                                                                                         |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| <b>Tested applications</b> | ICC, WB                                                                                                                                 |
| <b>Dilution range</b>      | WB=1:500-2000, ICC/IF=1:50-200                                                                                                          |
| <b>Antibody Type</b>       | Primary Antibody                                                                                                                        |
| <b>Clone Number</b>        | 2A3                                                                                                                                     |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>                | For research use only                                                                                                                   |



Sample: Lane 1: Muscle (Mouse) Lysate at 40 ug, Lane 2: Muscle (Rat) Lysate at 40 ug, Primary: Anti-Fbx32 (orb704278) at 1/1000 dilution, Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution, Predicted band size: 42 kD, Observed band size: 42 kD.

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)