

Product Datasheet

Goat Human IgG+IgM+IgA - Texas Red Antibody (orb705723)

Catalog Number	orb705723
Category	Antibodies
Description	The Goat Human IgG+IgM+IgA - Texas Red Antibody is suitable for ELISA, FC, IF, IHC. It is a Polyclonal, Texas Red conjugated antibody which raised against Human IgG+IgM+IgA. Purification: The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Clonality	Polyclonal
Species/Host	Goat
Conjugation	Texas Red
Reactivity	Human
Form/Appearance	1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Sodium Azide
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Immunogen	Human IgG+IgM+IgA
Tested applications	ELISA, FC, IF, IHC
Dilution range	IHC-P: 1:100-1:1000, IF: 1:100-1:1000, FC: 1:1000-1:4000
Antibody Type	Secondary Antibody
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)