

Product Datasheet

Rabbit Anti-Sheep IgG (H&L) Antibody (Cy5.5) (orb705876)

Catalog Number	orb705876
Category	Antibodies
Description	The Rabbit Sheep IgG (H&L) - Cy5.5 Antibody is suitable for ELISA, FC, IF, IHC. It is a Polyclonal, Cy5.5 conjugated antibody which raised against Sheep IgG. Purification: The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Clonality	Polyclonal
Species/Host	Rabbit
Conjugation	Cy5.5
Reactivity	Sheep
Form/Appearance	Liquid
Buffer/Preservatives	0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Sodium Azide
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Immunogen	Sheep IgG
Tested applications	ELISA, FC, IF, IHC
Dilution range	IHC-P: 1:100-1:1000, IF: 1:100-1:1000, FC: 1:1000-1:4000
Specificity	By immunoelectrophoresis and ELISA this antibody reacts specifically with Sheep IgG. No antibody was detected against non immunoglobulin serum proteins.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Antibody Type	Secondary Antibody
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)