

Product Datasheet

Goat Anti-Rabbit IgG (H&L) - Cy7 (orb705900)

Catalog Number	orb705900
Description	Goat polyclonal Secondary antibody against Goat Anti-Rabbit IgG (H&L) conjugated to Cy7
Species/Host	Goat
Reactivity	Rabbit
Conjugation	Cy7
Tested Applications	ELISA, FC, IF, IH
Immunogen	Rabbit IgG
Preservatives	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Sodium Azide
Form/Appearance	1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Antibody Type	Secondary Antibody
Source	Goat
Dilution Range	IHC-P: 1:100-1:1000, IF: 1:100-1:1000, FC: 1:1000-1:4000
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)