


Product Datasheet

Sucrose Phosphorylase Microplate Assay Kit (orb707336)

Catalog Number	orb707336
Description	Detection and Quantification of Sucrose Phosphorylase Activity.
Reactivity	All
Range	4 µmol/L - 400 µmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Tissue extracts, Cell lysate, Cell culture media, Other biological fluids, Platform: Microplate reader 340 nm
Sample Types	Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
CAS Number	9074-06-0
Expiration Date	Please enquire.



Representative Standard Curve

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)