


Product Datasheet

Sialic Acid Microplate Assay Kit (orb707350)

Catalog Number	orb707350
Description	Detection and Quantification of Sialic Acid Content.
Reactivity	All
Range	10 µmol/L - 1000 µmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Serum, Plasma, Cell culture media, Saliva, Milk, Tissue extracts, Cell lysate, Other biological fluids, Platform: Microplate reader 570 nm
Sample Types	Serum, Plasma, Cell culture media, Saliva, Milk, Tissue extracts, Cell lysate, Other biological fluids
Expiration Date	Please enquire.



Representative Standard Curve

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)