

Product Datasheet

Alpha Dystroglycan Rabbit Polyclonal Antibody (FITC) (orb7233)


Catalog Number	orb7233
Description	Alpha Dystroglycan Rabbit Polyclonal Antibody (FITC)
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	FITC
Tested Applications	IF
Immunogen	KLH conjugated synthetic peptide derived from human DAG1 (451-550/895aa)
Target	DAG1
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Antibody Type	Primary Antibody
Dilution Range	IF=1:200-800
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+441223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)



Paraformaldehyde-fixed, paraffin embedded (Mouse brain),
Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for
15 min, Block endogenous peroxidase by 3% hydrogen peroxide
for 20 minutes, Blocking buffer (normal goat serum) at 37°C for
30 min, Incubation: Anti-Alpha Dystroglycan Antibody,
conjugated (orb7233) 1:400, 1.5 hours at 37°C, DAPI (5 ug/ml,
blue) was used to stain the cell nuclei.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)