

Product Datasheet

E2F4 Antibody Blocking Peptide (orb72627)

Catalog Number	orb72627
Category	Proteins
Description	This product is E2F4 Antibody Blocking Peptide. It targets Transcription factor E2F4. Purification is performed using HPLC.
Target	Transcription factor E2F4
Form/Appearance	Lyophilized
Purification	HPLC
Storage	Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)