

## Product Datasheet

### Arr2 (Phospho-Ser366) Antibody (orb734357)

|                        |                                                                                           |
|------------------------|-------------------------------------------------------------------------------------------|
| <b>Description</b>     | Rabbit polyclonal antibody to Arr2 (phospho-Ser366)                                       |
| <b>Species/Host</b>    | Rabbit                                                                                    |
| <b>Reactivity</b>      | Insect, Plant                                                                             |
| <b>Conjugation</b>     | Unconjugated                                                                              |
| <b>Storage</b>         | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use. |
| <b>Note</b>            | For research use only                                                                     |
| <b>Clonality</b>       | Polyclonal                                                                                |
| <b>Antibody Type</b>   | Primary Antibody                                                                          |
| <b>Expiration Date</b> | 12 months from date of receipt.                                                           |

---

#### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

#### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)