

Product Datasheet

Rabbit Donkey IgG (H&L) Biotin Antibody (orb734546)

Catalog Number	orb734546
Category	Antibodies
Description	Rabbit Anti-Donkey IgG (H&L) Biotin
Clonality	Polyclonal
Species/Host	Rabbit
Conjugation	Biotin
Reactivity	Donkey
Antibody Type	Secondary Antibody
Storage	Store at -20°C/1 year
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)