

Product Datasheet

Mouse TNF Alpha/Tumor necrosis factor Quick ELISA Kit (orb738503)

Catalog Number	orb738503
Category	Assays and Kits
Description	Mouse TNF Alpha/Tumor necrosis factor Quick ELISA Kit
Target	Tumor necrosis factor
Reactivity	Mouse
Immunogen	Expression system for standard: E.coli; Immunogen sequence: L80-L235
UniProt ID	P06804
Tested applications	ELISA
Specificity	Natural and recombinant Mouse Tnf
Cross Reactivity	There is no detectable cross-reactivity with other relevant proteins.
Range	15.6 pg/ml - 1,000 pg/ml (cell culture supernatants), 7.8 pg/ml - 500 pg/ml (mouse serum, plasma)
Sensitivity	1 pg/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA)
Storage	Store at 4°C for 6 months. (Ships with gel ice, can store for up to 3 days in room temperature. Refrigerate upon receipt.)
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.


7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Mouse TNF alpha ELISA Kit

Mouse TNF alpha PicoKine ELISA Kit Standard Curve.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)