

Product Datasheet

DONKEY ANTI GOAT IgG (H/L):RPE (RAT/MOUSE ADSORBED) (orb746922)

Catalog Number	orb746922
Category	Antibodies
Description	Donkey polyclonal antibody against IgG conjugated to RPE
Target	IgG
Clonality	Polyclonal
Species/Host	Donkey
Isotype	Polyclonal IgG
Conjugation	RPE
Reactivity	Goat
Form/Appearance	Purified IgG conjugated to R. Phycoerythrin (RPE) - liquid
Concentration	0.5 mg/ml
Buffer/Preservatives	Purified IgG conjugated to R. Phycoerythrin (RPE) - liquid
Immunogen	Goat IgG.
Tested applications	FC
Storage	+4°C DO NOT FREEZE
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+441223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)