

Product Datasheet

GOAT ANTI RABBIT IgG (orb747119)

Catalog Number	orb747119
Category	Antibodies
Description	Goat polyclonal antibody to IgG
Target	IgG
Clonality	Polyclonal
Species/Host	Goat
Isotype	Polyclonal IgG
Conjugation	Unconjugated
Reactivity	Rabbit
Form/Appearance	Purified IgG - liquid
Concentration	IgG concentration 1.0 mg/ml
Buffer/Preservatives	Purified IgG - liquid
Purity	Purified
Immunogen	Purified rabbit IgG from normal rabbit serum
UniProt ID	P01870
Tested applications	DOT, ELISA
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)