

Product Datasheet

IgE Antibody (orb747864)

Catalog Number	orb747864
Category	Antibodies
Description	Mouse monoclonal antibody to IgE
Target	IgE
Clonality	Monoclonal
Species/Host	Mouse
Isotype	IgG1
Conjugation	Unconjugated
Reactivity	Animal
Form/Appearance	Purified IgG - liquid
Concentration	IgG concentration 1.0 mg/ml
Buffer/Preservatives	Preservative: 0.09% Sodium Azide (NaN ₃); Buffer: Phosphate buffered saline
Purity	Purified
Immunogen	Equine IgE
Tested applications	ELISA, IH-Fr, WB
Specificity	IgE
Antibody Type	Primary Antibody
Clone Number	3H10

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Storage


Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note


For research use only

Expiration Date

12 months from date of receipt.



Line graph illustrates about the Ag-Ab reactions using different concentrations of antigen and serial dilutions of IgG antibody.



Line graph illustrates about the Ag-Ab reactions using different concentrations of antigen and serial dilutions of IgG antibody.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558