

Product Datasheet

GOAT ANTI MOUSE IgG/A/M (H/L) (orb748907)

Catalog Number	orb748907
Category	Antibodies
Description	Goat polyclonal antibody to IgG IgA IgM
Target	IgG IgA IgM
Clonality	Polyclonal
Species/Host	Goat
Isotype	Polyclonal IgG
Conjugation	Unconjugated
Reactivity	Mouse
Form/Appearance	Purified IgG - liquid
Concentration	IgG concentration 1mg/ml
Buffer/Preservatives	Purified IgG - liquid
Purity	Purified
Immunogen	Purified mouse IgG, IgA and IgM.
UniProt ID	P01844, P01865, P01873, P01834, P01863, P01869, P03987, P01868, P01843, P01878, P01872, P01867, P01845, P01864
Tested applications	DOT, ELISA, FC
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Expiration Date 12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)