

Product Datasheet

GOAT ANTI MOUSE IgG/A/M (H/L):HRP (orb748908)

Catalog Number	orb748908
Category	Antibodies
Description	Goat polyclonal antibody against IgG IgA IgM conjugated to HRP
Target	IgG IgA IgM
Clonality	Polyclonal
Species/Host	Goat
Isotype	Polyclonal IgG
Conjugation	HRP
Reactivity	Mouse
Form/Appearance	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Buffer/Preservatives	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Immunogen	Purified mouse IgG, IgA and IgM.
UniProt ID	P01844, P01865, P01873, P01834, P01863, P01869, P03987, P01868, P01843, P01878, P01872, P01867, P01845, P01864
Tested applications	ELISA, WB
Storage	-20°C only (ship +4°C)
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)