

Product Datasheet

7-band Marker Protein (orb753179)

Catalog Number	orb753179
Description	8.5, 13.5, 18, 28, 40, 57 & 70 kDa Protein Marker
Preservatives	ProSpec 7-protein marker solution (0.5ug/ul/protein) contains 12mM Tris-HCl pH-6.8, 5% glycerol, 0.4% SDS 2.88mM 2-mercaptoethanol & 0.02% bromophenol blue.
Form/Appearance	Sterile Filtered blue solution.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles
Note	For research use only
Application notes	Recombinant & Natural Proteins
Source	E.Coli
MW	Use 5ul per lane. There is no need to heat before use. Each Protein conc. is 0.5ug/ul.
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)