

Product Datasheet

CoV-2 S1 (319-541) Protein (orb753211)

Description	The HEK293 derived recombinant protein contains the Coronavirus 2019 Spike Glycoprotein S1 Receptor Binding Domain [RBD], Wuhan-Hu-1 strain, amino acids 319-541 fused to His tag at C-terminal.
Preservatives	CoV-S1 RBD protein is 1mg/ml contains 20mM sodium phosphate (NaPP), 300mM NaCl, pH 7.2
Form/Appearance	Sterile Filtered clear solution.
Storage	Stability: CoV-2 S1 RBD Protein is shipped on ice packs. Upon arrival, Store at - 20°C
Note	For research use only
Application notes	Viral Antigens- Sars Proteins
Purity	Protein is > 90% pure as determined SDS-PAGE.
Source	HEK293 Cells
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)