

Product Datasheet

Human UBE2J2 protein (orb755556)

Catalog Number	orb755556
Category	Proteins
Description	Recombinant human UBE2J2 protein
Tag	N-terminal His Tag or N-terminal His-IF2DI Tag, determined during production process
Reactivity	Human
Buffer/Preservatives	Lyophilized from a 0.2 µm filtered solution in 10 mM Hepes, 500 mM NaCl with 5% trehalose, pH 7.4
Purity	Greater than 95% as determined by SDS-PAGE
UniProt ID	Q8N2K1
MW	11 kDa
Tested applications	ELISA, WB
Application notes	Centrifuge the vial at 10,000 rpm for 1 minute, reconstitute at 200 µg/ml in sterile distilled water by gentle pipetting 2-3 times, don't vortex
Source	E.Coli
Expression Region	124-224
Storage	-20°C for 12 months as lyophilized; 2-8°C for 1 month under sterile conditions after reconstitution
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558