

## Product Datasheet

### Human USP33 protein (orb755629)

<b>Catalog Number</b>	orb755629
<b>Category</b>	Proteins
<b>Description</b>	Recombinant human USP33 protein
<b>Tag</b>	N-terminal His Tag or N-terminal His-IF2DI Tag, determined during production process
<b>Reactivity</b>	Human
<b>Buffer/Preservatives</b>	Lyophilized from a 0.2 µm filtered solution in 10 mM HEPES, 500 mM NaCl with 5% trehalose, pH 7.4
<b>Purity</b>	Greater than 95% as determined by SDS-PAGE
<b>UniProt ID</b>	<b>Q8TEY7</b>
<b>MW</b>	25.5 kDa
<b>Tested applications</b>	ELISA, WB
<b>Application notes</b>	Centrifuge the vial at 10,000 rpm for 1 minute, reconstitute at 200 µg/ml in sterile distilled water by gentle pipetting 2-3 times, don't vortex
<b>Source</b>	E.Coli
<b>Expression Region</b>	71-303
<b>Storage</b>	-20°C for 12 months as lyophilized; 2-8°C for 1 month under sterile conditions after reconstitution
<b>Note</b>	For research use only
<b>Expiration Date</b>	6 months from date of receipt.

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)