

Product Datasheet

Human WFDC12 protein (orb755680)

Catalog Number	orb755680
Category	Proteins
Description	Recombinant human WFDC12 protein
Tag	N-terminal His Tag or N-terminal His-IF2DI Tag, determined during production process
Reactivity	Human
Buffer/Preservatives	Lyophilized from a 0.2 µm filtered solution in 10 mM HEPES, 500 mM NaCl with 5% trehalose, pH 7.4
Purity	Greater than 95% as determined by SDS-PAGE
UniProt ID	Q8WWY7
MW	9.6 kDa
Tested applications	ELISA, WB
Application notes	Centrifuge the vial at 10,000 rpm for 1 minute, reconstitute at 200 µg/ml in sterile distilled water by gentle pipetting 2-3 times, don't vortex
Source	E.Coli
Expression Region	24-111
Storage	-20°C for 12 months as lyophilized; 2-8°C for 1 month under sterile conditions after reconstitution
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)