

Product Datasheet

Human UBE2Z protein (orb756066)

| | |
|-----------------------------|--|
| Catalog Number | orb756066 |
| Category | Proteins |
| Description | Recombinant human UBE2Z protein |
| Tag | N-terminal His Tag or N-terminal His-IF2DI Tag, determined during production process |
| Reactivity | Human |
| Buffer/Preservatives | Lyophilized from a 0.2 µm filtered solution in 10 mM Hepes, 500 mM NaCl with 5% trehalose, pH 7.4 |
| Purity | Greater than 95% as determined by SDS-PAGE |
| UniProt ID | Q9H832 |
| MW | 27 kDa |
| Tested applications | ELISA, WB |
| Application notes | Centrifuge the vial at 10,000 rpm for 1 minute, reconstitute at 200 µg/ml in sterile distilled water by gentle pipetting 2-3 times, don't vortex |
| Source | E.Coli |
| Expression Region | 1-246 |
| Storage | -20°C for 12 months as lyophilized; 2-8°C for 1 month under sterile conditions after reconstitution |
| Note | For research use only |
| Expiration Date | 6 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)