

Product Datasheet

Human Kallistatin/Serpina4 Quick ELISA Kit (orb762975)

Catalog Number	orb762975
Category	Assays and Kits
Description	Human Kallistatin/Serpina4 Quick ELISA Kit
Target	Kallistatin
Reactivity	Human
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q21-P427
UniProt ID	P29622
Tested applications	ELISA
Specificity	Natural and recombinant Human SERPINA4
Cross Reactivity	There is no detectable cross-reactivity with other relevant proteins.
Standard	10 ng/tube
Range	156 pg/ml - 10,000 pg/ml
Sensitivity	10 pg/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA)
Storage	Store at 4°C for 6 months. (Ships with gel ice, can store for up to 3 days in room temperature. Refrigerate upon receipt.)
Note	For research use only

Biorbyt Ltd.


7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558



Human Kallistatin/Serpina4 PicoKine ELISA Kit standard curve.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)