

Product Datasheet

Human ANGPTL4/Angiopoietin Like 4/ARP4 Quick ELISA Kit (orb763002)

Catalog Number	orb763002
Category	Assays and Kits
Description	Human ANGPTL4/Angiopoietin Like 4/ARP4 Quick ELISA Kit
Target	Angiopoietin-related protein 4
Reactivity	Human
Immunogen	Expression system for standard: NS0; Immunogen sequence: L165-S406
UniProt ID	Q9BY76
Tested applications	ELISA
Specificity	Natural and recombinant Human ANGPTL4
Cross Reactivity	There is no detectable cross-reactivity with other relevant proteins.
Range	625 pg/ml - 40,000 pg/ml
Sensitivity	300 pg/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA)
Storage	Store at 4°C for 6 months. (Ships with gel ice, can store for up to 3 days in room temperature. Refrigerate upon receipt.)
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.


7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)



Human ANGPTL4 PicoKine ELISA Kit standard curve.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)