

Product Datasheet

HA-Tag Mouse mAb Antibody (orb763373)

Catalog Number	orb763373
Category	Antibodies
Description	Mouse monoclonal antibody to HA-Tag
Clonality	Monoclonal
Isotype	IgG1
Conjugation	Unconjugated
Concentration	1mg/mL
Buffer/Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Immunogen	Synthetic Peptide
Tested applications	IF, IP, WB
Dilution range	WB 1:5,000 IP 1:200 IF 1:1,000
Source	Mouse
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.


7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)


Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States


Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)



IP antibody use: 5ug HA Mouse IgG1 per ml Lysate, WB 1:5000
1: untransfected 293 cell lysate 2: transfected 293 cell lysate with HA-tag fusion protein 3: IP (untransfected 293 + anti-HA mAb + Protein G agarose) 4: IP (transfected 293 + normal Mouse IgG + Protein G ag).



IF analysis of 293 cells transfected with a HA-tag protein, using ELKbio anti-HATag Mouse mAb at a 1:2000 dilution (blue DAPI, red anti-HA).



0.5ug HA fusion protein + Primary antibody dilution at 1:10,000.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558