

Product Datasheet

Histone H2A.X(Acetyl Lys5) Rabbit pAb Antibody (orb763680)

Catalog Number	orb763680
Category	Antibodies
Description	Rabbit polyclonal antibody to Histone H2A.X (Acetyl Lys5)
Clonality	Polyclonal
Isotype	Antiserum
Conjugation	Unconjugated
Reactivity	Human, Mouse, Rat
Concentration	1mg/mL
Buffer/Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Immunogen	Synthetic Peptide
UniProt ID	P16104
MW	15kDa
Tested applications	IHC, WB
Dilution range	WB 1:1,000-2,000 IHC 1:200-500
Specificity	Human
Source	Rabbit

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Storage


Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note


For research use only

Expiration Date

12 months from date of receipt.



Western blot analysis of extracts from HeLa cells, untreated (-) or treated with TSA (1µM, 18 hr; +), using Acetyl-Histone H2A.X (Lys5) Rabbit pAb (upper) orb763680 1:2, 000 and β-actin Mouse mAb (lower) orb763396 1:5, 000.



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma using Histone H2A.X (Acetyl Lys5) (E116) Rabbit pAb diluted at 1:500.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558