

Product Datasheet

HIF-1 α Rabbit pAb Antibody (orb763844)

| | |
|-----------------------------|---|
| Catalog Number | orb763844 |
| Category | Antibodies |
| Description | Rabbit polyclonal antibody to HIF-1 α |
| Clonality | Polyclonal |
| Isotype | IgG |
| Conjugation | Unconjugated |
| Reactivity | Human, Mouse, Rat |
| Concentration | 1mg/mL |
| Buffer/Preservatives | PBS with 0.02% sodium azide and 50% glycerol pH 7.4. |
| Immunogen | Recombinant Protein |
| UniProt ID | Q16665 |
| MW | 120kDa |
| Tested applications | IHC, WB |
| Dilution range | WB 1:500-1,000 IHC 1:100-200 |
| Specificity | Human |
| Source | Rabbit |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)


Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States


Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Note For research use only

Expiration Date 12 months from date of receipt.



Immunohistochemical analysis of paraffin-embedded Human Colon Carcinoma Tissue using HIF-1 α (orb763844) Rabbit pAb diluted at 1:200.



Western blot analysis of 1) Mouse Liver Tissue Lysate, 2) Rat Liver Tissue Lysate, 3) HepG2 Cell Lysate using HIF-1 α (orb763844) Rabbit pAb diluted at 1:1000.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558