

## Product Datasheet

### PCBP1 Rabbit Polyclonal Antibody (orb77861)


|                            |                                                                                                                                         |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog Number</b>      | orb77861                                                                                                                                |
| <b>Description</b>         | Rabbit polyclonal antibody to PCBP1                                                                                                     |
| <b>Species/Host</b>        | Rabbit                                                                                                                                  |
| <b>Reactivity</b>          | Human                                                                                                                                   |
| <b>Conjugation</b>         | Unconjugated                                                                                                                            |
| <b>Tested Applications</b> | ELISA, IHC, WB                                                                                                                          |
| <b>Immunogen</b>           | Recombinant fragment corresponding to 64-325aa of human PCBP1.                                                                          |
| <b>Target</b>              | PCBP1                                                                                                                                   |
| <b>Form/Appearance</b>     | Lyophilized power with 1% rAlbumin and 0.02% NaN <sub>3</sub> .                                                                         |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>                | For research use only                                                                                                                   |
| <b>Application notes</b>   | Reconstitute in 100μL of sterile H <sub>2</sub> O with 50% glycerol.                                                                    |
| <b>Isotype</b>             | IgG                                                                                                                                     |
| <b>Clonality</b>           | Polyclonal                                                                                                                              |
| <b>MW</b>                  | 42KD                                                                                                                                    |
| <b>Dilution Range</b>      | WB: 1:200-1:1,000                                                                                                                       |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                         |

**Biorbyt Ltd.**


7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

**Biorbyt LLC.**

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)



Western blot analysis of Hela and K-562 cell lysate using PCBP1 antibody



Immunohistochemical analysis of formalin-fixed paraffin embedded fetal thymus showing cytoplasmic and nucleus staining using PCBP1 antibody

**Biorbyt Ltd.**

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

**Biorbyt LLC.**

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)