

## Product Datasheet

### ABCA7 Antibody (orb783919)

|                             |   |
|-----------------------------|---|
| <b>Catalog Number</b>       | orb783919   |
| <b>Category</b>             | Antibodies  |
| <b>Description</b>          | This ABCA7 Antibody is an unconjugated polyclonal product. It targets ABCA7 using a synthesized peptide derived from Escherichia coli O157:H7 aaeR as the immunogen. This antibody is suitable for ELISA, WB. Purification: Antigen Affinity. |
| <b>Target</b>               | ABCA7   |
| <b>Clonality</b>            | Polyclonal  |
| <b>Species/Host</b>         | Rabbit  |
| <b>Isotype</b>              | IgG   |
| <b>Conjugation</b>          | Unconjugated  |
| <b>Reactivity</b>           | A. thaliana   |
| <b>Form/Appearance</b>      | Liquid  |
| <b>Buffer/Preservatives</b> | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol   |
| <b>Purification</b>         | Antigen Affinity  |
| <b>Immunogen</b>            | Recombinant Arabidopsis thaliana ABCA7 protein  |
| <b>UniProt ID</b>           | <b>Q9STT5</b>   |
| <b>Tested applications</b>  | ELISA, WB   |

#### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

#### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Application notes**

1. Biorbyt guarantees the antibody purity > 95% confirmed by SDS-PAGE. 2. Biorbyt guarantees the antibody titer > 1: 64,000 confirmed by ELISA. 3. Biorbyt guarantees to offer antigen-specific Western Blot validation data. 4. Biorbyt guarantees 100% risk-free: No qualified antibody, no charge.

**Storage**

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

**Note**

For research use only

**Expiration Date**

12 months from date of receipt.

---

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)